

2015

# Cryo Crystallography

*Practical Course in X-Ray Diffraction: Module 3 + 4 (Protein Cryo Crystallography)*

In this module, the students will learn cryo sample preparation and cryo crystallographic data collection. This will involve harvesting a protein crystal from a crystal-tray using a CryoLoop, equilibrating the crystal in cryo protectant, flash cooling it in a cryogen like gaseous and/or liquid nitrogen, mounting the crystal in a goniometer and collecting x-ray diffraction data at cryogenic temperature of 100°K.



# Cryo Crystallography

---

## *Practical Course in X-Ray Diffraction: Module 3+4 (Protein Cryo Crystallography)*

### **Practical Course: When, What, Where?**

#### **Facility X-Ray Crystallography Facility**

Date	Wednesday, August 19, 2015
Event	Practical Course in X-Ray Diffraction Module 3+4
Venue	Kasha Lab Room Numbers 410 and 412

### **Module Objectives: Learn to perform cryo crystallographic experiment**

#### **Target In-coming graduate students**

Main Objective	The main objective is to learn how to harvest a protein crystal from a tray, flash cool it in a cryogen, mount the crystal in a goniometer and collect x-ray diffraction data
Other Objectives	Learning about cryogenics, cryogenic safety, flash cooling using liquid and gaseous cryogenics
	Learning about cryo protection, cryo buffers, and crystal handling
	Mounting a crystal in x-ray diffraction machine, aligning and collecting x-ray diffraction data

### **Resources: Book and web links**

#### **Subject Resource**

Book	Rhodes, G. (2006). <i>Crystallography made crystal clear</i> . Academic Press, New York.
Book	Blow, D. (2002). <i>Outline of Crystallography for Biologists</i> . Oxford University Press.
Web	UCLA Michael Sawaya's Notes on Cryo Crystallography: <a href="http://people.mbi.ucla.edu/sawaya/m230d/Data/data.html">http://people.mbi.ucla.edu/sawaya/m230d/Data/data.html</a>
Web	Elsbeth Garman Oxford England Lab's Help Tools: <a href="http://www.biop.ox.ac.uk/www/garman/lab_tools.html">http://www.biop.ox.ac.uk/www/garman/lab_tools.html</a>

### **Contact: X-Ray Facility Director**

#### **Instructor Thayumanasamy Somasundaram**

Office	414 Kasha Laboratory, 91 Chieftan Way, Florida State University, Tallahassee, FL 32306-4380
Phone	(850) 644-6448
E-mail	<a href="mailto:tsomasundaram@fsu.edu">tsomasundaram@fsu.edu</a>
URL	<a href="http://biophysics.fsu.edu/soma/projects/core-facilities-workshop/">http://biophysics.fsu.edu/soma/projects/core-facilities-workshop/</a>

## Table of Contents

Practical Course: When, What, Where? .....	1
Module Objectives: Learn to perform cryo crystallographic experiment.....	1
Resources: Book and web links .....	1
Contact: X-Ray Facility Director .....	1
Materials: List of chemicals and items required .....	3
How to do cryo crystallographic experiment.....	3
Cryo Crystallography of a protein crystal.....	4
§A. Cryo buffer preparation .....	4
§B. Harvesting a single crystal using a CryoLoop .....	5
§C. Cryo-buffer Treatment.....	5
§D. Flash-cooling in bulk liquid nitrogen .....	6
§E. Flash-cooling in gaseous nitrogen stream.....	7
§F. Aligning the Crystal .....	7

Version: 20150714 | 2015 IMB Boot camp | Copyright © 2012-15 Thayumanasamy Somasundaram

Version: 20150714; Original version: 20120628

Author: Thayumanasamy Somasundaram

Institute of Molecular Biophysics | Florida State University | Tallahassee, FL 32306-4380

Copyright © 2012-15 Thayumanasamy Somasundaram

# Cryo Crystallography

## Materials: List of chemicals and items required

### Chemicals Required:

Ethylene Glycol (Liquid)	Poly Ethylene Glycols (Solid & Liquid)	Sucrose (Solid)
Lysozyme Crystal	Ammonium Sulfate (Solid)	Sodium Chloride (Solid)

### Other Items Required:


Cryo Storage Dewar (10 L)	Cryo Foam Dewar (Short)	Cryo Foam Dewar (Tall)
CryoLoop®	CryoRod®	CrystalCap®
CrystalCap Magnetic®	Cryo Tongs°	CrystalWand®
Cryo Vial	Vial Clamp	Cover Slide
Forceps	Bottle Caps	

### Items for Crystal Set-up and diffraction:

Goniometer	Goniometer Magnet	Goniometer Wrench
Cryo Stream Blocker	Pipette Tips	Gilson® Pipettes
Optical Microscope	X-Ray Equipment	Crystal Scoring Sheet

## How to do cryo crystallographic experiment

Flash cooling and cryo crystallographic data collection of lysozyme crystal



The students will learn how to harvest a single-crystal (in this module the crystal of an enzyme lysozyme) grown using a technique called *hanging drop vapor-diffusion crystallization*. In this technique the crystals will grow in a drop of buffer hanging under the VDX plate's cover slip. A CryoLoop mounted on a CrystalCap will be used to pick up the crystal. The crystal will then be equilibrated in a cryo buffer. Then using either bulk liquid nitrogen or gaseous nitrogen maintained at 100°K (-173.15°C), the lysozyme crystal in the CryoLoop will be flash-cooled in milliseconds. The students will then mount the flash cooled crystal on goniometer attached to x-ray diffraction equipment and center the crystal. The students will then start x-ray diffraction data collection.

Cryo crystallography of protein crystal involves several steps and each step will be explained in detail next.